# RS Master Model Grey / Beige

## Gebrauchsinformation

## Deutsch

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## English

### Description:

**RS Master Model Grey** and **RS Master Model Beige** are light-curing resins for the additive manufacturing of high-precision dental models. The materials have been optimized for Rapidshape® DLP® D Series systems with a wavelength of 385 nm for layer thicknesses of 50 µm.

### Composition:

Composed of acrylates and initiators.  
**Please always observe the recommendations of the material safety data sheet!**

### Indication:

- Additive manufacturing of dental models, in particular prosthetic master models.

### Warning notes:

Avoid breathing vapours. Avoid contact with the unpolymerized liquid - wear protective gloves, protective clothing and eye protection. If on skin: wash with plenty of water. If skin irritation occurs, get medical advice. Uncured material is harmful to environment. Avoid release to the environment.

Always observe the applicable country-specific occupational health and safety regulations as well as the legal requirements when handling hazardous substances.

### Storage instructions / Shelf Life:

Protect uncured material from light, it may polymerize spontaneously. Store between 4°C and 25°C. Do not use after the expiration date. For shelf life, see label imprint.


### Disposal instructions:

Dispose of contents/container in accordance with local regulations.

Delivery Units:	REF
<b>D10+ Capsule</b> 10 x single use reservoir 10 x single use building platform 10 x RS Master Model Grey 60 g	RS001512
<b>D20+ Cartridge</b> 1 x reservoir 1 x RS Master Model Grey 200 g	RS001523
<b>RS Master Model Grey 1000 g</b>	RS001529
<b>D10+ Capsule</b> 10 x single use reservoir 10 x single use building platform 10 x RS Master Model Beige 60 g	RS001511
<b>D20+ Cartridge</b> 1 x reservoir 1 x RS Master Model Beige 200 g	RS001522
<b>RS Master Model Beige 1000 g</b>	RS001528

The product was developed for use in the dentistry and must be used according to the instructions for use. The Manufacturer does not accept any liability for damage caused by improper use of system or material. Furthermore, the user must independently check the material for its suitability and applicability for the intended purposes before using it, especially if this deviates from the process described in the instructions for use, or if this process is not described in the instructions for use.

## Application:

**Further processing information:**  
Parts made of **RS Master Model Grey** and **RS Master Model Beige** are subsequently processed and finished according to the normal dental technical processes.

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