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Instructions for use English

Description:
RS Model Ortho is a light-curing resin for the additive manufacturing of high-precision dental models. The material has been optimized for Rapidshape® DLP® D Series systems with a wavelength of 385 nm for layer thicknesses of 100 µm.

Composition:
 Composed of acrylates and initiators.
Please always observe the recommendations of the material safety data sheet!

Indication:

- Additive manufacturing of dental models, in particular for thermoforming models of orthodontic appliances such as aligners etc.

Warning notes:
 Avoid breathing vapours. Avoid contact with the unpolymerized liquid - wear protective gloves, protective clothing and eye protection. If on skin: wash with plenty of water. If skin irritation occurs, get medical advice. Uncured material is harmful to environment. Avoid release to the environment. Always observe the applicable country-specific occupational health and safety regulations as well as the legal requirements when handling hazardous substances.

Storage instructions / Shelf Life:
 Protect uncured material from light, it may polymerize spontaneously. Store between 4°C and 25°C. Do not use after the expiration date. For shelf life, see label imprint.

Disposal instructions:
 Dispose of contents/container in accordance with local regulations.


Delivery Units:	REF
D10+ Capsule 10 x single use reservoir 10 x single use building platform 10 x RS Model Ortho 60 g	RS001514
D20+ Cartridge 1 x reservoir 1 x RS Model Ortho 200 g	RS001712
RS Model Ortho 1000 g	RS001721

The product was developed for use in the dentistry and must be used according to the instructions for use. The Manufacturer does not accept any liability for damage caused by improper use of system or material. Furthermore, the user must independently check the material for its suitability and applicability for the intended purposes before using it, especially if this deviates from the process described in the instructions for use, or if this process is not described in the instructions for use.

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Application:

Further processing information:
 Parts made of **RS Model Ortho** are subsequently processed and finished according to the normal dental technical processes.

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