

MSE- Maxillary Skeletal Expander PRESCRIPTION



PLEASE PRINT

Account # C0 _____ **PO #** _____

B I L L I N G

PRACTICE TYPE: _____
(i.e., ortho, GP, pedo, prosth, oral surgeon, commer. lab)

DOCTOR: _____

ADDRESS: _____
(Specify if ship to address is different)

CITY: _____ STATE: _____ ZIP: _____

PHONE: (_____) _____

FAX: (_____) _____

EMAIL: _____

PATIENT: _____ AGE: _____

PLEASE PRINT

DATE SHIPPED: _____

DATE DUE: _____

1 day before appointment

LAB USE ONLY Incoming # cases _____

Customer Used: GLO Acct 2 Day On Call

Portal Upload - No Frt (00) Cust Acct - No Frt (00)

Disinfected 0 1 2 3 4 5 6 7 8 9

Rcvd: _____

B# _____ Via: _____

Shipment Date _____ Planned Shipment Date _____

(QC): _____ (LPD): _____

Estimated Delivery Date _____ Promised Delivery Date _____

ND _____ (Rec): _____

NO BITE / MDL - B / C Source: _____

Campaign: _____

Align ID# _____ Dig ID# _____

A D D R E S S

Please Provide: Boxes Labels

Rx: _____ Qty: _____
(specify appliance type)

IMPORTANT! Always retain models and bite until appliance is seated. Damage to models may occur during fabrication, please mark Rx if duplication (additional fee) of model(s) is required.

Maxillary Skeletal Expander:

Developed by Dr. Won Moon, is bi-cortically anchored to the palate with four TADs. Capable of dramatic skeletal expansion and nasal airway increase.

Standard MSE Features:

- Lab will determine the optimal MSE size
- Anterior TADs in line with mesial aspects of first molars
- Bands on first molars
- Screw arms anchored to first molar bands only
- Gap between screw and palate is <2mm
- 11mm TADs

Important:

- Doctor prefers an alternative A-P position of the MSE screw. Draw the line on model for the anterior TADs or use the line through mesial aspects of first molars for measurement reference.
- Please have a technical advisor contact me.

First MSE order - Includes the Starter Kit

Starter Kit Includes:

- (1) MSE Expansion Screw
- (4) Fixation TADs - 11mm
- (1) Ratchet Wrench Driver
- (1) Mini Driver Handle
- (1) Blade for Mini Driver/Handpiece
- (2) Spanner Wrenches (Activation)
- (2) Safety Leashes

Routine order - Includes the Refill Kit

Refill Kit Includes:

- (1) MSE Expansion Screw
- (4) Fixation TADs - 11mm
- (2) Spanner Wrenches (Activation)
- (2) Safety Leashes

MSE Expander Only

Kit Includes:

- (1) MSE Expansion Screw
- (2) Spanner Wrenches (Activation)
- (2) Safety Leashes

Bands: (Indicate band placement on illustration)


- Lab to provide Bands
- Lab to provide Bands with Buccal Tubes
- Lab to adapt doctors Bands

Protraction:

- Add Protraction Hooks

Additional Components:

- Mini Driver (Handle only) 281-052
- Mini Driver (Blade only) 281-030
- Ratchet Wrench Driver 281-053
- Long Spanner Wrench 281-050
- Short Spanner Wrench 281-051
- 11mm TAD (MSE-2) 281-025 Qty: _____
- 13mm TAD (MSE-2) 281-026 Qty: _____



License #: _____

Dr. Signature: _____

Master Rx on File # _____

Special Instructions: _____