

AM Aligner® Description


The AM Aligner a morning exercise tool to return the patient's bite back to its usual position after each use of oral appliance therapy. It is made of a unique thermoplastic material that is heated and fit to take a registration of the patient's bite before starting oral appliance therapy.

The AM Aligner can be used after use of any oral appliance on the market.



Key features:

- Tab allows for easier initial fitting and daily use
- ThermAcryl material is comfortable, resilient flexible, and thinner than the classic AM Aligner
- Shorter heating time and fitting time compared to the classic AM Aligner
- Anterior teeth guide helps patients bite down properly



How to fit

Fitting the AM Aligner takes less than 10 minutes. The new white color turns clear when ready to fit. Follow the instructions included in the packaging or [watch this fitting video](#).



Why use the AM Aligner®

All mandibular advancement devices can move the jaw forward permanently over time if the patient doesn't ensure that natural occlusion is reestablished each day.

The AM Aligner is to be used for 5-10 minutes each morning after use of an oral appliance to exercise the jaw back to its usual position.

Resources

[Link to AM Aligner fitting video.](#)