

Biostar® Scan with LCD Display



013-001

SCHEU
Dental Technology

Bar code scanner automatically enters the correct processing information for your material

Reaches working temperature in only **1** second.

Precision fabricate these appliances using positive pressure thermal-forming:

- Hawley Retainers
- Retainers with Bite Plane
- Expansion Appliances
- Bonded RPEs
- Invisible Retainers
- Custom Impression Trays
- Positioners
- Bleaching & Fluoride Trays
- TMJ Splints
- Temporaries for Crown & Bridge
- Indirect Bonding Trays
- Model Duplication
- Mouthguards
- Nightguards

Ask about our extensive line of high quality, break-resistant materials

Features:

- Advanced thermostatic infrared heating technology saves time and energy
- Uses up to 87psi positive pressure
- Extra deep model cup for embedding mounted models
- 3.2" color LCD display

Biostar Technical Information:

- Voltage: 110/120V - 50/60Hz
- Power: Maximum 750W
- Working Pressure: Max 6 bar/87psi
- W x D x H: 18" x 9.75" x 10.25"
- Weight: 35 lbs/14.5 kg



Great Lakes

ORTHODONTICS, LTD.

An Employee Owned Company

800.828.7626 • Fax: 716.871.0550 • greatlakesortho.com