



# keyprint<sup>®</sup>

in alliance with **LOCTITE**




## KeyModel Ultra™


KeyModel Ultra™ is a material designed for 3D printing of dental and orthodontic models.



**ULTRA-FAST TECHNOLOGY**



**FLAWLESS SURFACE DETAIL**



**THERMOFORMING RELEASE AGENT**



**IDEAL FOR THERMOFORMING ALIGNERS**

<p>● Ivory 1kg.....#285-431</p> <p>4kg.....#285-432</p>	<p>● Light Grey 1kg.....#285-429</p> <p>4kg.....#285-430</p>
---	--

PROPERTY	PROCEDURE	KEYMODEL ULTRA
Maximum Tensile Strength	ASTM D638	>50 MPa
Tensile Modulus	ASTM D638	> 1700 MPa
Elongation at Break	ASTM D638	5%
Maximum Flexural Strength	ASTM D790	>70 MPa
Flexural Modulus	ASTM D790	>1940 MPa
Viscosity @ 25C (cP)	ASTM D7867	500-600



**keystone**  
Industries



**Great Lakes**  
DENTAL TECHNOLOGIES

To order, please call 800.828.7626 or visit [GreatLakesDentalTech.com](http://GreatLakesDentalTech.com)