


KeyOrthoModel[®]


KeyOrthoModel is designed for high-speed and accurate prints. Perfect for thermoforming aligners and orthodontic devices.



ACCURATE



LOW SHRINKAGE



HIGH-SPEED

#285-427 (1kg bottle)

PROPERTY	PROCEDURE	KEYORTHOMODEL
Flexural Strength	ASTM D790	79.6
Flexural Modulus	ASTM D790	2293
Tensile Strength	ASTM D638	51
Tensile Young's Modulus	ASTM D638	1850
Elongation at Break (%)	ASTM D638	5.60%
Shore Hardness	Scale D	86

FEATURES

- OPTIMIZED TO PRINT @ 100 AND 200µm LAYER THICKNESS FOR HIGH THROUGHPUT
- PRINTS IN 1/2 THE TIME OF KEYMODEL AT SAME LAYER THICKNESS
- PAIRS WITH PRO-FORM[®]
- PERFECT FOR MAKING ALIGNERS
- WITHSTANDS THERMOFORMING TEMPERATURES



keystone[®]
Industries



Great Lakes
DENTAL TECHNOLOGIES

To order, please call 800.828.7626 or visit GreatLakesDentalTech.com