



# KeySplint Hard<sup>®</sup>

The newest member of the KeySplint family, KeySplint Hard is designed for rigid dental splints. KeySplint Hard provides therapy for bruxism and cases where tooth immobilization is needed.

## FEATURES

- BIOCOMPATIBLE MEDICAL DEVICE
- UNITED STATES FDA 510K-CLEARED
- HEALTH CANADA APPROVED
- CE MARKED (CLASS IIA)
- HARD FOR RIGID SPLINTS



**BIOCOMPATIBLE**



**DURABLE**



**EASY TO POLISH**



**EASILY CLEANED**



**ABRASION RESISTANT**

#285-421 (1kg bottle)

\*Currently available from multiple dealers. Please contact us

PROPERTY	STANDARD/METHOD	KEYSPLINT HARD
Flexural Strength	ISO 20795-2	60-65 MPa
Flexural Modulus	ISO 20795-2	1510-1600 MPa
Water Sorption	ISO 20795-2	18 µg/mm <sup>3</sup>
Elongation at Break	ASTM D638	9%
Biocompatibility	ISO 10993	Pass



**Great Lakes**  
DENTAL TECHNOLOGIES

To order, please call 800.828.7626 or visit [GreatLakesDentalTech.com](http://GreatLakesDentalTech.com)