

PLEASE PRINT

Account # LO _____ **PO #** _____

B I L L I N G PRACTICE TYPE: _____
(i.e., ortho, GP, pedo, prosth, oral surgeon, commer. lab)

DOCTOR: _____

ADDRESS: _____
(Specify if ship to address is different)

A D D R E S S CITY: _____ STATE: _____ ZIP: _____

PHONE: (____) _____

FAX: (____) _____

EMAIL: _____

PATIENT: _____ AGE: _____

PLEASE PRINT

DATE SHIPPED: _____

DATE DUE: _____
1 day before appointment

LAB USE ONLY	Incoming # cases	1	2	3	4	5+
Customer Used:	<input type="checkbox"/> GLO Acct	<input type="checkbox"/> 2 Day On Call				
<input type="checkbox"/> Portal Upload - No Frt (99)	<input type="checkbox"/> Cust Acct - No Frt (99)					
<input type="checkbox"/> Disinfected	0	1	2	3	4	5 6 7 8 9
Rcvd:	_____					
B#	_____	Via:	_____			
QC:	_____	LPD/Shp:	_____			
Needs DD Call	Rec:		_____			
NO BITE / MDL - B / C	Source:		_____			
Align ID#	_____	Dig ID#	_____			

Please Provide: Boxes Labels
 Rx: _____ Qty: _____
 (specify appliance type)
 Appliance Protection Program **(NOT AVAILABLE)**

**IMPORTANT! Always retain models and bite until appliance is seated.
 Damage to models may occur during fabrication, please mark Rx if duplication (additional fee) of model(s) is required.**

Medley Gold™ Sleep Apnea Appliance:

The Medley Gold appliance is 3 appliances in one. It allows the clinician to customize the same appliance with three different advancement mechanisms determined by patient needs: Rigid Nylon links, Elastomeric straps or Herbst* (Rod Sleeve). The Herbst version is PDAC approved.

Bite Registration- 5mm minimal vertical incisal clearance and 60% mandibular protrusion recommended. The laboratory may need to open the bite depending on curve of spee or in deep bite cases.

Medley Basic Appliance Platform Options:

Features the advancement location for Rigid Nylon Link or Elastomeric strap attachments. Includes elastic hooks, (1) package of elastics and (1) wrench.

Select preferred appliance

- Rigid Nylon Link - Includes sizes 17-22mm (Standard)
- Elastomeric - Includes sizes 16-21mm (Blue)

Medley Complete Appliance Platform (Rod Sleeve) PDAC Version

Features the advancement location for Telescopic Herbst arms (Rod Sleeve). Medicare Approved

Also includes: Rigid Nylon Link and Elastomeric strap attachments, (1) Rigid Nylon Link set (17-22mm), (1) Elastomeric strap set (16-21mm) - Blue, elastic hooks, (1) package of elastics and (1) wrench.

OPTIONS: Additional Fee

- Medley Configuration Kit** Includes: (1) set Rigid Nylon Links, (4) link screws, (1) set Elastomeric straps sizes (16-21mm) Blue & Yellow, (4) Elastomeric screws and (4) screw wrenches.

MATERIAL OPTIONS:

- Durasoft 2.5mm (Standard)
- Durasoft 1.8mm
- Biocryl with Thermacryl
- Hard Acrylic
- Clasp for Retention- Available only with Hard Acrylic material option:
 - Ball Clasps
 - Other _____

Recommended for use with Medley

- AM Aligner (10/pkg) 255-041

Master Rx on File # _____

Special Instructions: _____

For further details on the Medley and other sleep appliances visit our website at greatlakesdentaltech.com

Auxiliaries & Other Sleep Appliance Options:


Visit our website for a full list of **EMA Elastomeric Straps** available in various lengths and strengths in bulk packaging to meet your patient's needs.

- MPowRx™- Non-Custom Snoring and Sleep Apnea Appliance 255-040
- The George Gauge™ Kit 056-050
 - 2mm Bite Forks Grey Long (25/pkg) 056-032
 - 3mm Bite Forks Blue (25/pkg) 056-026
 - 5mm Bite Forks White Long (25/pkg) 056-033
- The George Gauge™ Asstd. Bite Fork Pack (12)-2mm, (12)-5mm 056-025
- my TAP™ 255-092
- AM Aligner (10/pkg) 255-041
- Mandibular Stabilizer _____ Red 155-028
- EZ Key (10/pkg) 140-020
- DentaSOAK® 235-023 (1 month Supply)
- Elastics (Bag of 50 Non-Latex) 833

Dental Midline Consideration: If no indication, we will follow bite registration.



Define Discrepancy of Dental Midlines:
 Upper to lower midline deviates to
 Pt. R or L by _____ mm



License #: _____

Dr. Signature: _____

*Herbst is a registered trademark of Dentaaurum, Inc.