


# PANTHERA | CLASSIC

YOUR PANTHERA CLASSIC NOW SHIPS WITH

## THE NEW 3.0 RODS


THE WIDEST RANGE OF RODS FOR THE WIDEST RANGE OF PATIENTS

- 20 PAIRS OF RODS PROVIDED IN 0.5 MM INCREMENTS WITH EVERY CASE
- MORE FREEDOM FOR TITRATION ADJUSTMENTS
- MORE TITRATION RANGE ACHIEVED WITH THE SAME SPLINT, NO REORDER NEEDED
- CUSTOM MADE RODS STILL AVAILABLE ON DEMAND



UNIQUE ANCHOR SYSTEM

LENGTH PRINTED ON THE ROD FOR EASY REFERENCE



ENGINEERED FOR OPTIMAL DEVICE RESISTANCE TO MANDIBULAR FORCES

TYPE 12 ORGANIC POLYAMIDE MATERIAL :

- Biocompatible, light & smooth
- Strong yet flexible & highly durable
- Does not elongate
- Identical material than the device: eliminates weaknesses at the connection points and reduces risk of fractures



**PANTHERA**  
SLEEP

800.828.7626 | GREATLAKESDENTALTECH.COM

Each 3.0 rod pack includes 2 x 10 pairs.