

TRIOS® 5 WIRELESS POD

The intraoral scanner designed for improved hygiene, optimal ergonomics, and effortless and faster precision scanning with ScanAssist.

—○ The 3Shape TRIOS® 5 IOS digital images are indicated for use in:

- Restorations
- Orthodontics
- Implantology
- Aid in diagnosis of caries



285-108

- ▶ The TRIOS® 5 Wireless Pod sets a new standard in patient protection and infection control. It is not only hygienic by design, but smaller, lighter, and designed to fit perfectly in every hand.
- ▶ It comes with the new **ScanAssist Engine** for perfect, faster scan results, and optional **TRIOS Share** for ultimate connectivity.
- ▶ **ScanAssist Engine** is intelligent alignment technology for building up a 3D model resulting in effortless and faster precision scans. Scan patients in any direction and create a scan path that works best. The new scan engine avoids stitching errors by reconstructing and building up the 3D model while scanning.
- ▶ **TRIOS Share** allows your entire practice to go digital. It is technology that lets you use your wireless TRIOS scanner on any PC in your practice. Now, your entire practice can be covered with just one scanner, so you can scan and plan treatment anywhere.

Note: TRIOS Share is FREE for all wireless scanners, however, a Splashtop subscription is required (sold separately).

TRIOS® 5 TECHNICAL SPECS AND INNOVATIONS

- 30% smaller and 20% lighter compared to TRIOS® 4
- Made to fit perfectly in your hand
- Introducing **ScanAssist**, providing easier scanning independent of scan strategy
- Improved battery technology, requiring only one battery per day
- Truly effortless, calibration-free scanner
- LED ring and haptic feedback for better user guidance
- Two-button design for improved software navigation and scan analysis
- Autoclavable scan tips providing a microbial barrier on the scan tip
- Body sleeves to significantly reduce risk of cross contamination
- Closed tip design to avoid cross contamination
- Scanner body was designed for cleaning and disinfection



More ▶

TRIOS® 5 WIRELESS POD



SCANNER SIZE

- Weight: 299g
- Light source: LED
- Output format: PLY, DCM, and STL
- Scanner dimensions: 26.6cm x 3.7cm x 3.8cm (LxWxH)
- Tip dimensions: 9.5cm x 2.26cm x 1.95cm (LxWxH)

WIRELESS RANGE

- Same room, up to 15 meters
- Scanning over WiFi possible with TRIOS Share

POD

- Diameter: 11cm
- Weight: 380g
- Place your TRIOS 5 on the pod in any orientation while not scanning.
- The pod is designed for easy cleaning and disinfection between patients for optimal hygiene.

BATTERY & USAGE

- (3) Rechargeable batteries
- Battery dimensions: 7.7cm x 2.17cm x 2.42cm (LxWxH)
- Battery duration: Instant heating and up to 66 minutes of continuous, non stop scanning per battery.
- Scanner sleep mode: Can stay in sleep mode with the battery inserted for more than 7 days.
- Charging time: 60 minutes to achieve 80% battery power. Less than 2 hours to achieve 100% battery power.

< INCLUDED in PACKAGE >

TRIOS 5 Scanner Wand
TRIOS 5 Scanner Pod
TRIOS 5 Battery Charger w/Power Cable
(5) TRIOS 5 Scanner Tips
TRIOS 5 Protection Tip
(100) TRIOS 5 Sleeves
(3) Batteries

Wireless Adapter
USB Dongle with Software
3Shape Subscription Dongle
Wall Mount Kit for Charger
3Shape Engagement Apps
Pod Installation Guide
Documentation

< Our Great Lakes Digital Specialists can review the benefits and options with you >
716.319.1260