

**PLEASE PRINT**

**Account # LO** \_\_\_\_\_ **PO #** \_\_\_\_\_

**B** PRACTICE TYPE: \_\_\_\_\_  
(i.e., ortho, GP, pedo, prosth, oral surgeon, commer. lab)

**I** DOCTOR: \_\_\_\_\_

**L** ADDRESS: \_\_\_\_\_  
(Specify if ship to address is different)

**A** CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

**D** PHONE: (\_\_\_\_) \_\_\_\_\_  
*Provide Country & City Code*

**R** FAX: (\_\_\_\_) \_\_\_\_\_

**E** EMAIL: \_\_\_\_\_

PATIENT: \_\_\_\_\_ AGE: \_\_\_\_\_

**PLEASE PRINT**

DATE SHIPPED: \_\_\_\_\_

DATE DUE: \_\_\_\_\_  
1 day before appointment

LAB USE ONLY	Incoming # cases	1	2	3	4	5+
Customer Used:	<input type="checkbox"/> GLO Acct	<input type="checkbox"/> 2 Day On Call				
	<input type="checkbox"/> Portal Upload - No Frt (99)	<input type="checkbox"/> Cust Acct - No Frt (99)				
<input type="checkbox"/> Disinfected		0	1	2	3	4
		5	6	7	8	9
Rcvd:	_____					
B#	_____	Via:	_____			
QC:	_____	LPD/Shp:	_____			
Needs DD Call		Rec:	_____			
NO BITE / MDL - B / C		Source:	_____			
		Location:	_____			
Align ID#	_____	Dig ID#	_____			

**Please Provide:**  Boxes  Labels

Rx: \_\_\_\_\_ Qty: \_\_\_\_\_  
(specify appliance type)

**Design Your Own**

**OR**

**Choose a New! eXceed Kit**

**eXceed Rx-** brackets are placed on the pretreatment models using a standard placement method. Bracket placement can be modified using the eXceed software and approved by the doctor.

**eXceed Tx-** brackets custom positioned using a simulated final setup and virtual arch wire. The setup is based upon clinical instructions from the customer. Bracket positions are selected to optimize this treatment and can be modified by the customer using the eXceed software. Patient photos and Xrays are required.

**Bracket Information:**

- Brackets Enclosed
- Bracket Inventory at the Lab
- Use GLO house brackets
  - EasyClip+
  - BioTru Classic
  - BioTru Ceramic
  - BioTru Sapphire
  - Add Molar Tubes

**Slot Size:**  .018 or  .022

**Prescription:**  Roth  MBT  Damon (EasyClip+ Only)

**Bracket:**  Upper  Lower  Both  
 5 - 5  6 - 6  7 - 7

**EasyClip+ Kit**  
eXceed Rx + EasyClip+ brackets + BioTru Contour Nitinol Archwires

**BioTru Classic Kit**  
eXceed Rx + BioTru Classic twin brackets + BioTru Nitinol Archwires

**BioTru Ceramic Kit**  
eXceed Rx + BioTru Ceramic brackets + BioTru Nitinol Archwires

**Aesthetic Express Kit**  
eXceed Rx + BioTru Ceramic brackets + Aesthetic Short Term Archwires

**Please send the following supplies:**

- Reliance Flow-Tain® Composite (065-040)
- Reliance Assure® Plus Bonding Resin (065-046)
- Reliance Liquid Etching Agent (065-004)
- Elastics \_\_\_\_\_

**Tray Fabrication Options:**

**Dual Tray** (hard outer shell, soft inner shell) for all orders.  1 Piece (Standard)  Midline Split  3 Pieces Per Arch

Master Rx on File # \_\_\_\_\_

Special Instructions: \_\_\_\_\_

License #: \_\_\_\_\_

Dr. Signature: \_\_\_\_\_