Code 253

Heating Time: **110 sec.**

DARK BLUE BIOCRYL
2.5mm/125mm Round
021-008

MINISTAR S ▪ SCANNER

BIOSTAR VII ▪ SCANNER

Heating Time
Easy-to-Scan
Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes
DENTAL TECHNOLOGIES

800.828.7626 ▪ greatlakesdentaltech.com

SMPL471Rev061121