Code 143

**Heating Time: 40 sec.**

**PURPLE MOUTHGUARD**
**2mm/125mm Round**
**021-033**



**MINISTAR S ▪ SCANNER**



**BIOSTAR VII ▪ SCANNER**

Heating Time
Easy-to-Scan
Barcode Strip

**NOTE:** Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes™
DENTAL TECHNOLOGIES

**800.828.7626 ▪ GLDT.com**

SMPL508Rev032526