Code 233

Heating Time: **90 sec.**

RED MOUTHGUARD
3mm/125mm Round
021-038

MINISTAR S ▪ SCANNER

BIOSTAR VII ▪ SCANNER

Heating Time
Easy-to-Scan
Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes
DENTAL TECHNOLOGIES

800.828.7626 ▪ greatlakesdentaltech.com

SMPL513Rev093021