Code 233

Heating Time: 90 sec.

YELLOW MOUTHGUARD
3mm/125mm Round
021-044

MINISTAR S ▪ SCANNER

BIOSTAR VI ▪ SCANNER

GREAT LAKES ORTHODONTICS

Heating Time Easy-to-Scan Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

800.828.7626 ▪ greatlakesortho.com

GREAT LAKES ORTHODONTICS

SMPL519Rev040417