Code 143

Heating Time: **40 sec.**

**ASSORTED MOUTHGUARD**
**PURPLE**
**2mm/125mm Round**
**021-045**



**MINISTAR S ▪ SCANNER**



**BIOSTAR VII ▪ SCANNER**

Heating Time
Easy-to-Scan
Barcode Strip

**NOTE:** Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes
DENTAL TECHNOLOGIES

**800.828.7626 ▪ greatlakesdentaltech.com**

SMPL524Rev012320

Code 183

Heating Time: **60 sec.**

**ASSORTED MOUTHGUARD**
**GREEN, RED, ORANGE,**
**BLUE, YELLOW, WHITE**
**2mm/125mm Round**
**021-045**



**MINISTAR S ▪ SCANNER**



**BIOSTAR VII ▪ SCANNER**

Heating Time
Easy-to-Scan
Barcode Strip

**NOTE:** Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes
DENTAL TECHNOLOGIES

**800.828.7626 ▪ greatlakesdentaltech.com**

SMPL525Rev012320

Code 123

**Heating Time: 30 sec.**

**ASSORTED MOUTHGUARD**
**BLACK**
**2mm/125mm Round**
**021-045**



**MINISTAR S ▪ SCANNER**



**BIOSTAR VII ▪ SCANNER**

Heating Time
Easy-to-Scan
Barcode Strip

**NOTE:** Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes
DENTAL TECHNOLOGIES

**800.828.7626 ▪ greatlakesdentaltech.com**

SMPL526Rev012320