Code 233

Heating Time: 90 sec.

ASSORTED MOUTHGUARD
GREEN, RED, ORANGE,
BLUE, YELLOW, WHITE
3mm/125mm Round
021-046

MINISTAR S ▪ SCANNER

BIOSTAR VII ▪ SCANNER

Heating Time
Easy-to-Scan
Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes™
DENTAL TECHNOLOGIES

800.828.7626 ▪ greatlakesdentaltech.com

SMPL528Rev030922

Code 203

Heating Time: 70 sec.

ASSORTED MOUTHGUARD
PURPLE
3mm/125mm Round
021-046

MINISTAR S ▪ SCANNER

BIOSTAR VII ▪ SCANNER

Heating Time
Easy-to-Scan
Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes™
DENTAL TECHNOLOGIES

800.828.7626 ▪ greatlakesdentaltech.com

SMPL527Rev030922

Code 143

Heating Time: 40 sec.

ASSORTED MOUTHGUARD
BLACK
3mm/125mm Round
021-046

MINISTAR S ▪ SCANNER

BIOSTAR VII ▪ SCANNER

Heating Time
Easy-to-Scan
Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes™
DENTAL TECHNOLOGIES

800.828.7626 ▪ greatlakesdentaltech.com

SMPL529Rev030922