Code 183

Heating Time: 60 sec.

LAMINATE MOUTHGUARD
PINK GLITTER
4mm/125mm Round
021-122

MINISTAR S ▪ SCANNER

BIOSTAR VII ▪ SCANNER

Heating Time
Easy-to-Scan
Barcode Strip

NOTE: Heating times listed serve as a general guideline only. Individual results may vary according to factors such as model size or room temperature. You may need to adjust heating times according to your individual situation.

Great Lakes™
DENTAL TECHNOLOGIES

800.828.7626 ▪ greatlakesdentaltech.com

SMPL688Rev091522